


COSIATEC and SIATEC Compress: Pattern Discovery by Geometric Compression




David Meredith dave@create.aau.dk

1. Representing music with point sets




• Piece of music represented as a set of points in pitch-time space

2. Maximal translatable patterns (MTPs)



• Maximal translatable pattern (MTP) for a given vector in a given dataset contains all points that can be translated by that vector to other points in the dataset

3. Translational equivalence classes (TECs)



- Translational equivalence class (TEC) of a pattern contains all translationally invariant occurrences of the pattern in a given dataset
- A translational equivalence class can be compactly expressed as a (pattern, vector set) pair, giving the points in one occurrence of its pattern along with the vectors that map that occurrence onto all the other occurrences
- The **covered set** of a TEC is the union of all the occurrences in the TEC (e.g., the coloured points in each of the diagrams above)

$\mathbf{T} \leftarrow \langle angle$

6. The COSIATEC algorithm


```
COSIATEC(D)
           P \leftarrow \text{Copy}(D)
           T^* \leftarrow \mathsf{nil}
           while P \neq \emptyset
                   T^* \leftarrow \text{GetBestTEC}(P, D)
                   \mathbf{T} \leftarrow \mathbf{T} \oplus \langle T^* \rangle
                   P \leftarrow P \setminus \text{COV}(T^*)
           return T
```

- COSIATEC (COmpression with SIATEC) partitions a dataset into the covered sets of a set of MTP TECs
- SIATEC is run on the dataset and the 'best' TEC is selected and added to the output (lines 5–6)
- The best TEC is the one that gives the best compression. Compactness is used as a tie-breaker.
- The best TEC's covered set is removed from the dataset (line 7)
- The process is repeated until there are no points left (i.e., $P = \emptyset$ in line 4)
- The output is a sequence of TECs, T, such that the covered sets of these TECs form a partition of the input dataset
- By representing each TEC as a (pattern, vector set) pair, the output encoding is shorter than the explicitly encoded input dataset

7. The SIATECCompress algorithm


- Computes the MTP TECs using SIATEC
- Sorts TECs into decreasing order of 'quality'
- Quality of a TEC determined by compression ratio, with compactness of pattern used as a tie breaker
- Descends list of sorted TECs, choosing TECs that cover sufficient number of uncovered points
- Encoding consists of a list of TECs, in decreasing order of quality, such that their covered sets cover the input dataset
- Covered sets of TECs in output may intersect (unlike COSIATEC)
- Encoding produced by SIATECCompress typically less compressed than that produced by COSIATEC

4. The SIA algorithm




- SIA (Structure Induction Algorithm) discovers all the MTPs in a set of n, k-dimensional points in $O(kn^2 \log_2 n)$ time and $O(kn^2)$ space in the worst case
- Finds vector from each point to every lexicographically later point and stores vector in a table with a pointer back to the origin point
- (vector, point) pairs sorted lexicographically, then segmented at points where vector changes
- Points in each resulting segment form an MTP

8. Raw, BB and Segment versions




- Three versions of each algorithm (COSIATEC and SIATECCompress) submitted to MIREX 2013 competition
- Raw versions: patterns found are the raw MTP occurrences (blue points in left figure)
- BB versions: patterns found are the contents of the bounding boxes of the MTP occurrences (blue points in middle figure)
- Segment versions: patterns found are the contents of the segments spanning the MTP occurrences (blue points in right figure)

5. The SIATEC algorithm



- SIATEC (SIA+TECs) discovers all MTP TECs in a kdimensional dataset of size n in $O(kn^3)$ time and $O(kn^2)$ space in the worst case
- Executes SIA, but fills in complete vector table: set of vectors by which an MTP is translatable is intersection of columns in vector table headed by points in the MTP

9. Results on MIREX 2013 Test Database



- DM10 (SIATECCompressSegment) generally performed best overall with respect to precision, recall and F1 measure.
- NF2 and DM10 not significantly different on establishment recall on a per-pattern basis and NF2 significantly faster than DM10.
- Full results available at http://www.music-ir.org/mirex/wiki/ 2013:Discovery_of_Repeated_Themes_&_Sections_Results